



ADVY

Partnering with you



EzDx™ | AdvDx™

Rapid Malaria Antigen Detection Test Kits



WHO PQ Approved*



CE Marked



NIMR Evaluated



The Global Fund*

Where Precision meets Performance.

EzDx™ | AdvDx™

Malaria Rapid Test

Features

- In-house developed highly accurate pLDH and HRP-II Monoclonal Antibodies.
- Highly sensitive and specific.
- Results in 20 minutes.
- Storage temperature from 2°C to 40°C.
- Excellent stability.


Highlights

- Easy to use
- 24 months shelf life at 2°C - 40°C.
- Sample volume 5 µL.
- Whole blood protocol.
- Clear background for interpretation.

EzDx™ MALARIA Pv/Pf Ag

Plasmodium falciparum (Pf) and Plasmodium vivax (Pv) antigen


Test Procedure



Product Performance

	Sensitivity	Specificity
Pf	97.77%	98.42%
Pv	98.33%	98.42%


Result Interpretation



EzDx™ MALARIA Pan/Pf Ag

Plasmodium falciparum (Pf) antigen and P. vivax, P. ovale, or P. malariae antigen


Test Procedure



Product Performance

	Sensitivity	Specificity
Pf	97.77%	97.89%
Pv	97.50%	97.89%


Result Interpretation



AdvDx™ MALARIA Pf Ag

Plasmodium falciparum (Pf) antigen


Test Procedure



Product Performance

	Sensitivity	Specificity
Pf	96.66%	98.42%

Result Interpretation



*applicable to AdvDx™ only

EzDx™

Malaria Rapid Test



AdvDx™

Malaria Rapid Test



PRODUCT ORDERING INFORMATION:

Product Name	Product Code	Description	Pack size
EzDx™ Malaria Pan/Pf Antigen	00-DKM-RK-MAL-EDX-001-001	Detection of Pan(pLDH) & Pf(HRP-II) Antigen	1 Test / Kit
	00-DKM-RK-MAL-EDX-001-010		10 Tests / Kit
	00-DKM-RK-MAL-EDX-001-025		25 Tests / Kit
	00-DKM-RK-MAL-EDX-001-050		50 Tests / Kit
EzDx™ Malaria Pv/Pf Antigen	00-DKM-RK-MAL-EDX-002-010	Detection of Pv(pLDH) & Pf(HRP-II) Antigen	10 Tests / Kit
	00-DKM-RK-MAL-EDX-002-025		25 Tests / Kit
	00-DKM-RK-MAL-EDX-002-050		50 Tests / Kit
AdvDx™ Malaria Pf Ag	00-DKM-RK-MAL-ADX-004-025	Detection of Pf(HRP-II) Antigen	25 Tests / Kit

ANTIGENS | ANTIBODIES | DIAGNOSTIC KITS

+91 865 7428 612
 info@advychemical.com
 www.advychemical.com



@advy
 @advy_dx
 @advydx
 @advychemical